WILLI GELLER
Creation VC
W O R K B O O K

WWW.CREATION-WILLIGELLER.COM

WILLI GELLER
Creation VC

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Composite is and can more than you expect!


Due to its high filler content, composite is a material that can bear comparison with dental ceramic with regard to durability and aesthetic performance; in some respects, for example the opalescence, it is even far superior. Why do so many dental technicians still shy away from this material? Why does composite live a shadow existence – absolutely unjustified?

*Experience in this workbook that composite is more and can more! Our authors passionately dealt with Creation VC – everybody in his or her own way – and revealed the significance, the benefit and the possibilities of this modern hybrid composite. On the following pages they show you what characterises Creation VC and how versatilely and aesthetically it can be used – from **material science** (Robert Arvai) to the **basic layering** (Nando Aeschlimann) to the **individual layering** (Uwe Gehringer) and from the **press technique** (Annette v. Hajmasy) to the **combined restoration** with Creopal Shells and Creation VC Gingiva (Michael Zangl).*

*What distinguishes Creation VC compared to other composites? On the one hand, it follows the market requirements, on the other hand the high standards we claim for all our products. Thus Creation VC impresses not only with its **outstanding aesthetics and functionality**, but also with its **versatile, flexible use**. Similar to the Creation porcelains, VC offers very opaque, colour intensive dentines and opaque dentines which can be combined with the various enamel, effect and transpa materials to achieve natural looking veneers even with thin layers. And as Creation VC is available in the variants Paste and Flow, it offers a huge playground of processing techniques.*


*In short: Creation VC offers a **great potential for dental technology** and more than complies with our brand promise – “simply brilliant”! Be inspired and excited by our authors and Creation VC with this workbook!*

Thomas Hiebel
Marketing & Sales Manager

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


88–105

Creopal Shell, Creation VC /
Creopal Shell, Creation VC
Michael Zangl

CREATION VC **PASTE** – MATERIALKENNWERTE / CREATION VC **PASTE** – MATERIAL PROPERTIES

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Creation VC pastes were developed according to the hybrid microfiller concept (MFR) and feature a reinforced structure of organic and inorganic fillers in a low-resin bonding matrix.

Creation VC offers high mechanical values and provides excellent physical properties. The results are very smooth surfaces, a good polishing behaviour, a high wear resistance, biocompatibility and an antagonist-friendly abrasion behaviour.

Vickers hardness (HV): 35

Flexural strength: 124.1 MPa

Silicate filler portion: 7.5 %–69.6 % (* this value varies for the different materials)*

CREATION VC **FLOW** – MATERIALKENNWERTE / CREATION VC **FLOW** – MATERIAL PROPERTIES

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Biegefestigkeit: 165 MPa

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Specially developed flowable composite with outstanding material properties covering a wide field of applications due to its high flowability.

The material composition is based on microfillers made of superfine glass fillers with a particle size of 200 nm showing exceptionally high resistance to high masticatory forces.

With patients suffering of bruxism, this composition prevents breaking and chipping of particles at high forces. An additional positive effect is the low refractive index resulting in a very natural translucency.

The latest silane system and the high content of microfillers of 69.9 % provide a high homogeneity without clusters.

Creation VC possesses high strength values and a high wear resistance due to the extremely high bond of the particles in the composite matrix.

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
1 x Creation VC Cervical Opaque (2,4 ml): O-CO

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DESCRIPTION OF MATERIALS

1 x Creation VC Base Opaque (2.4ml): O-BO

Paste opaque with good flow properties and excellent light-curing qualities. Dyed in a bright, pale yellow colour and also polymerises in areas with undercuts. Serves as bonder for the opaques O-A1–O-D4.

1 x Creation VC Cervical Opaque (2.4ml): O-CO

The Cervical Opaque is applied after the VC Base Opaque along the gingival line with a width of 1 mm. It covers metallic margins and can be used to characterise e.g. fissure areas of posterior teeth.

16 x Creation VC Opaque (2.4ml): O-A1–O-D4

The paste opaques in 16 A–D basic tooth shades can be applied easily and possess an excellent opacity.

5 x Creation VC Opaque Dentine (2.9ml): OD-32, -37, -41, -43, -44

Opaque dentines can be used instead of the shades of the standard range or be mixed with them to achieve diverse effects out of the depth of the crown.

- **OD-32 havanna:** warm shade of yellow for individual shading in the dentine, incisal or fissure area.
- **OD-37 curry:** the yellow shade is well suitable for the interproximal area to improve the light reflection and to achieve a warm shade effect.
- **OD-41 orange:** the yellow-orange shade is very well suitable for the application to the oral surfaces of the restoration. It provides liveliness without being too glaring or too dark.
- **OD-43 ivory:** light ivory. Especially suitable for areas where high brightness values have to be achieved. These materials can be used e.g. in the middle third of a crown to increase the brightness.
- **OD-44 cuba:** dark and intensive shade for application in the cervical area and to reproduce exposed root parts.

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DESCRIPTION OF MATERIALS

16 x Creation VC Dentine (2.9 ml): D-A1 – D-D4

The development of the 16 dentine pastes was adjusted to the Vita® shade guide. Due to the opacity and the chroma of the Creation VC dentines, it is possible to apply only very thin layer thicknesses. This is a special advantage when only minimal space is available. The recommended layer thickness of the dentine is between 0.5 and 0.7 mm and can be modified according to experience.

4 x Creation VC Enamel (2.9 ml): E-57 – E-60

The scope of delivery includes four standard enamels: E-57, E-58, E-59 and E-60. Shades from white to grey can be achieved with these materials and they can be used in different concentrations for diverse effects.

E-57 (very high brightness), E-58 (high brightness), E-59 (medium brightness), E-60 (low brightness). The concept of the colour brightness and the experience of the dental technician offer a wide range for the application of these enamels depending on the desired effect. E-57 can be used for example to brighten the dental ridge or cusp tips. E-60 on D-A2 for example produces a greyer incisal effect.

2 x Creation VC Pearl Enamel (2.9 ml): PS-1, PS-2

The pearl incisal materials are suitable to imitate demineralisations or other white discolorations. PS-1 and PS-2 can be used to accentuate white areas or for cusp tips.

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3 x Creation VC Effect Enamel (2,9ml): SI-01, SI-02, SI-03

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DESCRIPTION OF MATERIALS

The SI-Effect Enamel and the TI-Transpa Enamel materials were developed to offer the dental technician even more possibilities to reproduce the enamel shades. These unique materials have opalescent characteristics and appear more vivid than standard enamel materials. They can be used pure or mixed with the basic enamels.

3 x Creation VC Effect Enamel (2.9ml): SI-01, SI-02, SI-03

The three SI-Effect Enamels are similar to the basic enamels (E-57–E-60) regarding opacity (capacity of light reflection). Each material has a particular yellow-orange shade intensifying from SI-01 to SI-03. The SI-Effect Enamels have an identical opalescence independent of the colouring. An increase of the colour saturation reduces the brightness. A light SI-01 at the dental ridge or the cusp tips results in a bright area appearing warm and vivid. In contrast, SI-03 can be used on the occlusal surfaces of elderly people to naturally reproduce the darker enamel area.

4 x Creation VC Transpa (2.9ml): TI-01, TI-02, TI-04, TI-05

TI-01 blue: *creates a gentle, translucent, blue effect (opacity) at the mesial and distal incisal edge of an anterior restoration – similar to natural teeth.*

TI-02 white: *accentuates cusp tips or dental ridges in the posterior area and results in a bright, matt appearance. This material can also be used to accentuate mesial and distal line angles of an anterior tooth. In principle, the use of a thin layer of TI-02 between the dentine and the enamel is recommended to achieve a finely nuanced effect.*

TI-04 yellow: *can be used almost unlimitedly for every restoration. When TI-04 is applied undiluted at the lingual area of the anterior, the material gives the crown a warm, translucent and vivid appearance.*

TI-05 grey: *creates a distinct depth effect and can be used for the reproduction of old or very transparent teeth.*

MASSEN BESCHREIBUNG

3 x Creation VC Neck Transpa (2,9ml): HT-52, HT-53, HT-54

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DESCRIPTION OF MATERIALS

3 x Creation VC Neck Transpa (2.9ml): HT-52, HT-53, HT-54

The Neck Transpa materials can be used in the cervical area as well as in all other areas of the restoration. The name "Neck Transpa" derives from the highly translucent colour which can often be seen at the tooth neck of natural teeth.

These materials can create interesting shade hues depending on the placement on the restoration. The following examples describe possible applications of the material which can be developed further in practice. The materials can be used pure as well as mixed with dentines or enamels.

HT-52 khaki: pure in the incisal area of an anterior restoration or as wash material in the occlusal area.

HT-53 sand: pure or mixed with enamel material especially suitable for the occlusal or incisal area.

HT-54 honey: to achieve darker shades, this colour can be mixed with enamel materials and be used in the incisal, cervical or occlusal area of the restoration.

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9 x Creation VC Flow Dentine, 2 x Flow Enamel, 2 x Flow Opaque Dentine (2,0ml):

FD-BW | FD-A1, FD-A2, FD-A3, FD-A3.5, FD-B1, FD-B2, FD-B3 | FD-HT | OD-A2, OD-A3, EN-L, EN-T.

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DESCRIPTION OF MATERIALS

Creation VC Flow: optimal handling and first-class aesthetics

Creation VC flow is a first-class composite with excellent stability, polishing characteristics and versatility which can be used for restorations of all kind.

Its unique filler technology results in a high filler content and a uniform distribution of the fillers. Creation VC Flow is best suitable for the creation of indirect composite restorations and it is the perfect completion to Creation VC Paste. In addition, the Creation VC Flow materials can perfectly be combined with the Creopal teeth and the Creopal Shell veneering shells. Thus they allow the realisation of modern therapy concepts such as the creation of implant-supported prostheses. Very complex restorations can consequently be created forecasted and efficiently. The results are durable tooth restorations with outstanding aesthetics.

The specially designed ergonomic tip with coating allows an easy application with controlled dosage of the flow materials. In the injection technique, they can be injected into a hollow mould like a flowable composite, but they can also be layered like a classical composite.

The Creation VC Flow materials can be mixed and combined with Creation VC Make up Neo for an individual adjustment and a better aesthetic result.

9 x Creation VC Flow Dentine, 2 x Flow Enamel, 2 x Flow Opaque Dentine (2.0ml):

FD-BW | FD-A1, FD-A2, FD-A3, FD-A3.5, FD-B1, FD-B2, FD-B3 | FD-HT | OD-A2, OD-A3, EN-L, EN-T.

Universal, light-curing, radiopaque and flowable composite. Creation VC Flow is a multi-indicative composite with controlled flow rate.

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Stains and modifiers for individual characterisation

Creation VC Make up Neo stains offer a large number of stains for the characterisation so that perfect aesthetics are achieved. Make up Neo is produced with the proven nano-filler technology and provides great abrasion resistance and thus long-lasting gloss. With the complete “two-in-one” stain system, individual shade accentuations can be added using the internal or external staining technique.

- Versatile and efficient
- High shade accuracy and fidelity due to first-class materials
- Easy handling: can be mixed among each other and are correctable

Fields of application:

- Mixed in for intensification or individualisation
- Depth effect with MUN-IL Illusion
- Tooth shaded colouration of the restorations

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12 x Creation VC Make up Neo (2,6 ml): MUN-1 – MUN-12, MUN-IL, MUN-G

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1 x Creation VC Make up Neo (2,6 ml): Translucent h.v. (high viscosity)

1 x Creation VC Make up Neo (2,6 ml): Translucent l.v. (low viscosity)

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Application steps for staining or characterisation:

Internal staining with MUN:

1. Optionally also mixing with Translucent h.v. or l.v. for individual adjustment of the translucency (increase) and consistency during layering.
2. Subsequent over-layering of the fixed shades.
3. Preparation for glazing resp. polishing.
4. Completion of the restoration.

3 x Creation VC Make up Neo (2.6 ml): MUN-A, -B, -C

The stains MUN-A to MUN-C are used to intensify the chroma in the dentine and tooth neck area.

12 x Creation VC Make up Neo (2.6 ml): MUN-1 – MUN-12, MUN-IL, MUN-G

The 12 MUN stains can be added in small quantities to modify the Creation VC Flow materials for an increase of the colour saturation. The stains MUN-IL (Illusion) and MUN-G (Gingiva) are suitable for special effects at crowns and gingiva.

1 x Creation VC Make up Neo (5 ml): MUN-CL

Colour neutral material to mix with MUN stains.

1 x Creation VC Make up Neo (5 ml): MUN-CL h.c. (high consistency)

Colour neutral material with same effect as Make up Neo MUN-CL, with increased viscosity.

1 x Creation VC Make up Neo (2.6 ml): Translucent h.v. (high viscosity)

1 x Creation VC Make up Neo (2.6 ml): Translucent l.v. (low viscosity)

Composites with different consistencies to mix (dilute) the MUN stains.

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Creation VC Gingiva is excellently suitable for the lifelike reproduction of the gingiva tissue.

The assortment is available in three consistencies:

- Opaque (Creation VC Gingiva Opaque),
- Paste (Creation VC Gingiva) and
- Flowable (Creation VC Gingiva Modifier and Translucent Gel).

With its opaques and the layering materials in two different consistencies, the light-curing composite enables manifold options for a lifelike reproduction of the oral mucosa. Veins can be imitated with fine fibres to be added into the veneering. The anatomic principles have to be obeyed for a reproduction close to nature. For example, the area of the keratinised gingiva has to be designed in light pink due to the normally reduced blood circulation in this part. In contrast, however, the mucogingival area has a strong blood circulation and is intensely veined.

4 x Creation VC Gingiva Opaque (2.4ml): O-G3, O-G4, O-G5, O-G6

Gingiva-shaded opaques to cover the framework structures prior to veneering with Creation VC Gingiva or Creation VC Modifier. The opaques have a high opacity and are easy to process.

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4 x Creation VC Gingiva (2.9ml): G3, G4, G5, G6

Composite in paste-like consistency. The shades were developed according to the requirements of natural, multi-cultural gingiva shades. Creation VC Gingiva can be used separately or in combination with Creation VC Gingiva Modifiers.

6 x Creation VC Gingiva Modifier (2.4ml): GM-1, GM-2, GM-3, GM-5, GM-6, GM-7

Creation VC Gingiva Modifier is a flowable composite in 6 shades. Due to its low viscosity and opacity, it can be mixed very well and be used for the individual characterisation. The shades vary over hues red and cover with GM-1 purple the complete scope of shades necessary to reproduce the gingiva.

1 x Creation VC Gingiva Translucent (2.4ml)

Flowable, translucent composite. This composite is especially suitable to be mixed with the Gingiva Modifiers to soften their shades.

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1 x Creation VC Separator

The Separator is an insulation plaster to resin. It is applied onto the plaster model prior to the creation of inlays and onlays.

1 x Creation VC Air Protect

This product is applied on the composite veneering prior to the final polymerisation to prevent the formation of an inhibition layer and to guarantee a complete polymerisation of the composite.

1 x Creation VC Diamond Polisher

The diamond-containing polishing paste was developed especially for Creation VC and accounts for a high shine and a brilliant appearance of the veneerings.

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The Composite Primer serves as bonding agent between the microfilled hybrid composite Creation VC and other resins like e.g. microfilled composite or acrylics and can be used for corrections at Creation VC composite veneerings or for repairs of acrylics.

1. Blast resin surface with 2–3 bar aluminium oxide (50–110µm).
2. Clean restoration with 70 % isopropyl alcohol in an ultrasonic unit at 70 °C or with dry, oil-free compressed air.
3. Apply one to two drops VC Composite Primer with a brush as a thin layer onto the surfaces to be bonded.
4. Cure the Creation VC Composite Primer with a light-curing unit (wave length: 470 nm) for one minute.
5. After light-curing, do not remove the oxygen inhibiting outer layer which serves as bonding layer for Creation VC veneering composites and stains.
6. Apply veneering composite and cure.

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The silane-containing bonding agent creates a chemical bond between metal and composite. The silanes of the Metal Primer react chemically with the Creation VC Base Opaque and thus ensure the safe bond between metal and the composite Creation VC.

1. Blast metal surface with 2–3 bar aluminium oxide (50–110µm).
2. Clean restoration with 70 % isopropyl alcohol in an ultrasonic unit at 70 °C or with dry, oil-free compressed air.
3. Apply Creation VC Metal Primer with a clean brush onto the retention surface.
4. Let dry at 60 °C. Suitable are an infrared lamp, a pre-heating furnace or a drying cabinet.
5. Apply a thin layer (wash) Creation VC Base Opaque onto the retention surface with a flat brush. The Base Opaque serves as bonder between the veneering framework and the subsequent layers. The further layering is carried out with opaques or layering materials.
6. Cure every layer with a light-curing unit.

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The Creation VC Ceramic Primer is a bonding agent for the conditioning of the adhesion surfaces of glass ceramic, zirconium, aluminium oxide and hybrid porcelains.

- 1. Blast bonding surface with 100µm aluminium oxide (1–3 bar according to the framework material).
- 2. Clean with 70 % isopropyl alcohol in an ultrasonic unit at 70 °C or with dry, oil-free compressed air.
- 3. Apply Creation VC Ceramic Primer in a thin layer with a brush and blow excess with compressed air.
- 4. Further processing with Creation VC Base Opaque as bonding agent and in the next step with the layering materials.

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*Composites are biocompatible, antagonist-friendly and can be repaired and extended. They are **versatile with controlled and comfortable handling** as well as comfortable for the patient. And the photo-optical characteristics of composites are closer to the natural tooth as those of porcelains. Nevertheless, they are hidden in the darkest drawer in most laboratories! Maybe composites are often marginalised because they are considered to be “difficult”, a necessary evil. But they are not difficult; you only have to know how they are composed and how you have to treat them correctly.*

*My former employment in the industry aroused my interest in material science. I understood that the material production is based on science and there is no alternative to reading books and studies. Later I became a member of the “Bio Emulation” group and then the great library with evidence based knowledge was at my disposal. In dental technology, there are maybe many ways leading to Rome, but only one is appropriate for the material involved. Going that way, many mistakes can be avoided from the very beginning. Only the exact know-how of the materials allows us to leave the path and try something new. So: **composites are not difficult**; they are reliable hybrid all-rounders if you treat them correctly:*

1. Correct surface treatment of the framework.
2. Accurate priming and bonding.
3. Correct light-curing.
4. No steam cleaning.

And why Creation VC? Because it is a high-performance material due to its composition. It can be combined and used in many ways, offers a range of brilliant colours, can be polished well and features high hardness. I am a passionate user of Creation CC and VC also fills a gap as it is perfectly adjusted to the colour range of the metal ceramic.

*We as dental technicians should engage more in composite and put it on top of our material list. The material is a **jack-of-all-trades** that not only upgrades our spectrum but may also lead to a professional exchange at eye level with the dentist – since he uses composite every day.*

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


Bowen RL (1956), *Use of Epoxy Resin in Restorative Materials*. J. Dent. Res (35:360-369).

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Material science of dental composites

Nowadays a multitude of the most different materials is used in dental, oral and craniomandibular sciences and there is no end in sight. The fields of application of the materials are manifold, thus every patient requiring artificial dentition will get in contact with dental materials. A material which has been indispensable in dental laboratories for many years is the dental composite. First publications date back until the year 1956 and were published by Dr. Rafael Bowen.

Bowen RL (1956), Use of Epoxy Resin in Restorative Materials. J. Dent. Res (35:360-369).

Composition of composites


Composites are materials composed of two or more different components and chemical phases. Different substances are embedded in an organic carrier material (matrix) and have a considerable influence on the characteristics. Such substances can be fillers, colour pigments, ultraviolet stabilisers or radiopaque additives. Compared to the self-curing and chemically curing composites from the pioneer era originally developed as an alternative to amalgam, today's materials are hybrid all-rounders used in a much wider field of application.

Creation VC is a composite at the latest state-of-the-art and fulfils all requirements of the processing in the dental laboratory.

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






The resin matrix is also called organic phase and is the carrier resp. base material of the composite.

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Functions of the organic phase:

- Bonds all ingredients of the composite in one compound
- Can be polymerised (responsible for the material formulation)
- Is the carrier material for pigments and fillers
- Controls the contraction

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Fillers are called the dispersed phase. They develop by homogeneously mixing filler materials like e.g. silicone dioxide with a liquid monomer and polymerisation with heat. In the next process step, the cured resin with the dispersed filler materials is crushed and ground to a powder. The filler particles have a size of 5 to 50 micrometers. The silanes developing at the surface form a very good bond with the organic phase. The filler portion of a modern dental composite is approx. 73 % (percentage by volume).

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Functions of the dispersed phase:

- Improvement of the mechanical characteristics of the composite (e.g. flexural strength, polishing properties)
- Increase of the wear resistance
- Reduction of the polymerisation shrinkage
- Reduction of the coefficient of thermal expansion

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Coupling phase

In the coupling phase, the metastable laboratory composite is polymerised into a stable state with the influence of light. During this process, the light initiators included in the matrix, e.g. camphor quinone or benzoin alkylether, are activated initiating the polymerisation.

The illustration shows how the organic matrix and the fillers bond to each other. The bond develops by mechanic connection (rough surface) and by chemical reaction (silane-carbon bond).

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Fig. 7: Coupling phase during the light polymerisation of composites.

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The key to success for restorations from composites is the making of good bond to the framework and the correct processing.

Surface treatment and surface preparation:

- Sandblasting of the framework with 110 µm at 3–4 bar
- Cleaning with isopropyl alcohol in an ultrasonic unit for 5 minutes (no steam cleaner)
- Addition of the silane primer (Creation VC Metal Primer) immediately after the evaporation of the alcohol. This primer can also be used with PEEK or PEKK
- Application of the bonder according to the instructions of the manufacturer (Creation VC Base Opaque)
- Correct light-curing of the composite (Creation VC Base Opaque)
- Correct surface treatment with milling cutters and polishers after curing
- Proper cleaning (no steam cleaner)

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Effect of Incorrect Primer Application on Dentin Permeability

Alexandre Cavalcheiro/Marcos A. Vargas/
Steven R. Armstrong/Deborah V. Dawson/
David G. Gratton

A review of adhesion science

Sally J. Marshall, Stephen C. Bayne, Robert
Baier, Antoni P. Tomsia, Grayson W. Marshall

The effect of different surface treatments on bond strength between leucite reinforced feldspathic ceramic and composite resin

Andrea Fabianelli, Sarah Pollington, Federica
Papacchini, Cecilia Goracci, Amerigo Cantoro,
Marco Ferrari, Richard van Noort

Effect of drying time of 3-methacryloxypropyltrimethoxysilane on the shear bond strength of a composite resin to silica-coated base/noble alloys

Mutlu Özcan, Jukka P. Matinlinna, Pekka K.
Vallittu, Marie-Charlotte Huysmans

Longevity of resin composite restorations

Shisei Kubo, Center for Clinical Education and
Training, Nagasaki University Hospital, 1-7-1
Sakamoto, Nagasaki 852-8588, Japan

Randomized Clinical Trial of Indirect Resin Composite and Ceramic Veneers: Up to 3-year Follow-up

Marco M.M. Gresnigt/Warner Kalk/Mutlu
Özcan

Properties of an indirect composite material polymerized with two different laboratory polymerizing systems

Dent Mater J. 2005;24(3):377-81. Satsukawa
H, Koizumi H, Tanoue N, Nemoto M, Ogino T,
Matsumura H.

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There is not the one and only standard layering as teeth are as different as humans. Every tooth is a little unique piece of art with its shape, shade and structure. In addition, it changes over the years. A young tooth is not only brighter resp. shows juvenile mamelons, it also shows a more distinct surface texture. An older tooth however was used thus showing a high gloss level and incisal abrasion.

It is necessary to incorporate as much information and characteristics as possible into a reconstruction so that it can develop from the crown towards the tooth. Thus the ambitious dental technician moves away from a basic layering to an individual solution. If the ratio dentine/enamel is correctly adjusted and the natural tooth characteristics are integrated, a good result can be achieved even with a simplified layering.

Creation VC behaves remarkably similar to nature with regard to light optics, therefore the fundamentals of a tooth can be reproduced with only few materials and parameters. An authentic result can thus be achieved with minimal effort.

On the following pages I describe **two age-specific basic layerings** (junior and senior tooth) following the basic tooth structure: dentine core and incisal are designed with a few resources analogue to the internal morphology of natural teeth.

For the **juvenile tooth**, I used only two materials and concentrated on the composition of the dentine mamelon. I used no special effects as the opalescence already develops with the enamel material. For the **senior tooth**, I added a stain internally to make it appear older from the depth. With my case documentation I particularly want to show how telescopic works can be created with a transparent silicone key and composite and how easy but still close to nature long-term temporaries can be composed with only a few light-dynamic materials.

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The shape and the size of teeth can be reproduced with silicone keys. The fully anatomical impression is therefore taken from the labial and palatal side. In addition, the silicone key serves as control of the space available when material is applied for the individual characterisation. After curing, this first silicone key of the palatal dental arch is insulated with a silicone insulation.

The second key from labial is produced with a transparent silicone as hard as possible from the cartridge. The desired tooth shape after the individualisation of the dentine core can thus be reconstructed without difficulties. The composite can be pre-polymerised through the silicone key due to its translucency.

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Even with a basic layering, the structure of the dentine core significantly influences the characteristic of a tooth, so the development should be executed very accurately. So take the school pencil again and sketch the

reduction of the dentine on the fully anatomical restoration. The dentine core determines the brightness and the chroma as well as the layout of the mamelons and the degree of opalescence.

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Dentine cut-back:

Consider the following notions: is the incisal material darker or lighter than the dentine core and where should the tooth show the lightest area? In this case, the enamel material (E-59) is darker than the dentine of tooth 11 and lighter than the dentine core of tooth 21. The more the incisal area of tooth 21 is reduced, the darker it will appear by applying the enamel material E-59.

The opposite is the case for tooth 21: brightness will arise with the enamel material where dentine is reduced. Tooth 21 shall have the highest brightness in the medium third, therefore only very little of the fully anatomical crown is reduced in this area. The tips of the mamelons are also reduced only very slightly in labial direction so that they are later visible.

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The finally reduced dentine cores are blasted with aluminium oxide 25–50µm and 1.5bar. Clean the restoration only with oil-free compressed air or isopropyl alcohol. In the next step, apply Creation VC Composite Primer with a microbrush in a thin layer onto the cleaned and dry surface. Blow off the

excess with compressed air. Then cure the VC Composite Primer with light. An inhibition layer develops guaranteeing the bond for the layering of further materials.

Abb. 12: Composite Primer für den Verbund.
Fig. 12: Composite Primer for bonding.

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Minimalism: only the darker dentine core is incisally characterised with Creation VC Make up Neo (MUN-4). The tooth will have an age-appropriate appearance with only little effort.

Technical procedure:

After the polymerisation of Creation VC Composite Primer, shake well the bottle with Creation VC Make up Neo. Give some drops into a mixing bowl. Apply Creation VC Make up Neo on tooth 21 incisal-centrally in the upper third onto the dentine surface using an instrument or a brush. Light-cure Make up Neo immediately after application.

Tip: Creation VC Make up Neo was used pure for the case described. It can be diluted with Creation VC Translucent h.v. and l.v. for the internal characterisation of for subtle hues.

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Fig. 16: Characterisation with MUN-4 (mandarine). Staining is not necessary for the juvenile tooth.

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
Fig. 17: Creation VC cartridge.

- Enamel E-59

The enamel material E-59 is very well suitable for the pressing-over of both dentine cores to achieve the desired appearance of the crowns. The material emits a certain warmth and possesses the optimal translucency. As the dentine core of the juvenile tooth was not individualised and the tooth shape was only slightly reduced, it can be assumed that the silicone keys can be positioned on the model without any problems. The individualisation in

terms of colour also creates no significant volume increase at the labial area. For pressing the Creation VC pastes should be heated slowly. Suitable are heat lamps, warm water or a blow-drier. The viscosity of the paste is reduced by warming, the risk of the inclusion of air bubbles is reduced and the excess material can be displaced more easily during the pressing.

Attention:

Creation VC must never be heated above 60 °C! The model, the two silicone keys and the reduced restorations can also be heated up prior to pressing resulting in a better flowability. For the pressing itself, the transparent silicone key is lined bubble-free with tempered E-59 at ambient light. Subsequently, the filled silicone key is placed slowly into its final position on the model. Polymerisation is

effected in a light-curing unit after the correct position was reached and the pressing excess is very little.

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Fig. 19: Visible excesses after pressing and pre-polymerisation.

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For the complete polymerisation of the veneerings, the matrix is removed and Creation VC Air Protect is applied onto the surfaces. This blocker prevents the forming of an inhibition layer on the polymerised restoration. Afterwards the final polymerisation takes place with the light-curing units GC LABOLIGHT LV-III or GC LABOLIGHT DUO* for 3 minutes. The Creation VC Air Protect gel is washed off with water. Final corrections can now be executed by grinding or polishing of the*

surfaces. Absolutely avoid an overheating of the restoration as the composition of the veneering composites and thus the physical properties change with temperatures of 60 °C and above.

** Customary light-curing units with comparable light power can also be used.*

Abb. 20: Air Protect-Gel.
Fig. 20: Air Protect gel.

Abb. 21: Restauration mit Air Protect zur Endaushärtung.
Fig. 21: Restoration with Air Protect for final curing.

Abb. 22: Creation Diamond Polisher.
Fig. 22: Creation Diamond Polisher.

Abb. 23: Finale Restaurationen (jugendlicher und gealterter Zahn).
Fig. 23: Final restorations (juvenile and aged tooth).

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*It is not courage but especially **clarification of and communication with the patient** to create individual aesthetics. If the dental technician can exchange views with the patient concerning wishes and possibilities, the two are on the right path for individual tooth aesthetics in more than 90 % of the cases.*

*As long as we design tooth defects like disharmonies, discolourations and enamel cracks or noticeable tooth shapes and positions in a way that they result in a **harmonic appearance**, the new restoration will mostly be incorporated invisibly. Nothing is the same in nature, it is always individual – and it is important to imitate nature with all its peculiarities.*

*Thus, patient-oriented tooth aesthetics are not possible for me and my workflow without personal patient contact. I have to know in advance if the patient likes an individual layering as much as I do. Only when I found that out and we **speak the same language** I can start with the creativity and realisation process and approach the final restoration with wax-ups, mock-ups and prototypes until everything fits.*

The aesthetic individual result of the following case could only be achieved by an intensive and repeated interaction with the patient. The basis for our communication about the decision-making was a prototype with set-up Creopal teeth. This was the only way for the patient to know what his new anterior teeth would look like and that individual aesthetics would suit him much better than a design too uniform and too regular.

*The **individual layering of the maxilla anteriors** – the effects are very well shown in the photos – can hardly be noticed in the patient's mouth within speaking distance: the new teeth look as natural as if they were natural! Creation VC helped me a lot as it offered me all possibilities for an individual layering technique and processability with its brilliant stains and effect materials, the high processing quality and the various consistencies – also resp. especially in limited spaces.*

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

Fig. 3: Application of five opaque layers in the posterior tooth area as preparation of the polychrome pressing over (1st layer: OB-O, 2nd to 4th: O-A2, 5th: O-B3 occlusal).

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Fig. 4: The posterior teeth pressed over three times (1st: OD-37, 2nd: D-A2, incisal application: TI-01, 3rd E-59) with flat occlusal surfaces matching the mandible not yet restored.

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Fig. 5: Rubbing in of the Metal Primer with a brush, let dry for some seconds.

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Fig. 6: First non-covering opaque layer – wafer-thin to guarantee curing (light-curing time: 1 min.).

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Fig. 9: Fifth, darker opaque application (approximal and at the canine teeth) to achieve maximum vitality from the depth.

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Fig. 10: Application of different opaque dentines to control opacity and warmth in the restoration.

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Fig. 13: The flow dentine is slightly more translucent than the paste dentine, thus a smooth transition from dentine to incisal can be achieved.

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Fig. 14: Application of the mamelons to control the incisal translucency (OD-37 + G5 in the mixing ratio 1:1).

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Fig. 15: Application of Flow Enamel EN-T incisal at the tooth neck with a transparent matrix, then creation of the ridges and overlaying of the body with incisal materials.

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Fig. 16: Sharpening of the blue ridges with a separating disc to lay out the enamel cracks at the transition enamel to dentine.

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Fig. 20: Completion of the incisal tooth shape with incisal materials and neck transpa material HT-52 (above E-59).

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Fig. 21: Reduction of layering flaws or bubbles (comparable to the practice after the first dentine firing when processing layering ceramic).

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Fig. 22: Incorporation of additional enamel cracks with a separating disc for a more complex appearance of the tooth structure.

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*Cooking depends on the ingredients, the selection of spices and the experience and craftsmanship of the cook – all that decides whether the dish tastes well. It is the same with us: we need a good plan and a material perfectly suiting the present case and offering a range of shades and effect materials for the **individual “flavour”** of our work. And of course we need to be experts in our handicraft and have to know how to process which ingredients and which nuances are possible. In short: we need a coherent recipe!*

With the following case I want to show how important it is to think in advance about the selection of techniques and materials and to exactly define shape, length and position of the teeth.


*When the patient is satisfied with the diagnostics, I use the **pressing technique with composite** to implement the final restoration. Using this technique I can absolutely maintain the shape developed together with the patient and confirmed by him. Forward-looking I can implement the promise I made in the trial run 1:1. Additionally I can compact the composite and do not have to worry about later bubbles and air inclusions.*

*The challenge of the present case was the need to combine abutment teeth with implants. Milled PEEK copings were glued into the tertiary framework as compensating elements to compensate possible tensions between the immovable implants and the easily movable abutment teeth. As this **triple construction** requires a lot of space I decided to use Creation VC with its very opaque, chromatic dentine. I wanted to try whether I can create an aesthetic veneering with depth effect with this micro-hybrid composite – despite limited spaces. After the dentine pressing and cut-back I did not think in units and glossaries but like a good cook in colours. I mixed nuances, placed effects, added gingiva materials and thus **individually refined my work** – exactly according to the patient's taste!*

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Due to medical reasons only one-part zirconium implants with their typical telescope-like structure could be used to increase the number of abutments.

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Fig. 2: Initial situation in the mandible: teeth 33 and 37 could be preserved.

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Fig. 3: The mandible was implanted at the tooth positions 34, 36, 43 and 46.

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Fig. 4: Initial situation in the maxilla: teeth 21 and 22 could be preserved.

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Fig. 5: The maxilla was implanted at the tooth positions 11, 13, 16 and 26.

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Procedure:

the restoration of the natural abutment teeth is created with telescopic primary crowns, the insertion direction is parallel to the telescopic structure of the implants. PEEK copings (0.2 to 0.3 mm) are produced as secondary structure.

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Fig. 8: Control of the placed PEEK primary crowns in the articulator.

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The non-precious metal frameworks are produced as reduced form of the fully anatomical resin bridges.

For the try-in, the framework is glued with the PEEK copings in the mouth (“passive fit”) and a new impression is taken to create the new working model.


Abb. 9: Tertiärstruktur der Oberkieferversorgung.
Fig. 9: Tertiary structure of the maxilla restoration.


Abb. 10: Tertiärstruktur der Unterkiefferversorgung.
Fig. 10: Tertiary structure of the mandible restoration.

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Fig. 11: The so-called dummies are created in resin with the available 3D data of the temporary restoration; the dummies are transferred to the new model situation. The final correction is effected in wax. The set-ups serve as basis for the creation of the composite veneerings.

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Fig. 12–14: Resin bridge invested with silicone as fully anatomical spacer. You save the wax-up on the framework as the dummy is simply substituted by the metal framework after the investment.

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Fig. 20 and 21: To protect the glued in PEEK friction parts, they are covered with silicone at the internal areas prior to blasting.

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
At first Creation VC Metal Primer is applied onto the blasted frameworks. In the next step the first opaque layer is applied very thinly as wash opaque. The very light, colour-neutral Creation VC Base Opaque is suitable and it is polymerised for 60 seconds.

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
Fig. 22 and 23: Blasted frameworks conditioned with Metal Primer and with Creation VC Base Opaque.

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78

The clear silicone counter is now filled with the respective dentine of Creation VC. The flask halves are put together and closed with the three fixation screws for the pressingover of the opaqued framework. Now the restoration is light-cured.

Abb. 26: Füllung des Klarsilikon-Konters mit Creation VC-Dentinpaste D-A2.
Fig. 26: Filling of the clear silicone counter with Creation VC dentine paste D-A2.

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After curing in the light-curing unit, the flask is opened and the restoration taken out. For a complete curing of the basal portions of the restoration, it is cured a second time in the light-curing unit. Afterwards the pressing excess is removed.

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Fig. 28: Opened flask after pressing over.

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Prior to the individualisation, the processed surfaces must be conditioned with the Creation VC Composite Primer. It is thinly applied with a brush and then polymerised.

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After the individualisation, the restorations are placed back into the flask and pressed over with Creation VC Enamel pastes.

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Fig. 31: Individual intermediate layering – "internal application".

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Fig. 32: Incisal pressing after individualisation.

- Transpa TI-01 blue, TI-04 yellow
- Neck Transpa high fluorescent HT-53 sand
- Pearl Enamel PS-1 blue
- Enamel E-60

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The gingival areas of the restoration are modelled in wax to be poured later with PMMA. Only the area of the anterior teeth is designed individually with Creation VC Gingiva.

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The basic parts of the gingiva are produced in PMMA. Prior to the individualisation with Creation VC Gingiva, a bonding layer has to be created with Creation VC Composite Primer.

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How far am I allowed to go? How individually can I finish? How much “naturalness” should be seen? Every patient has a certain idea how he or she wants to look like – and it is my daily challenge to match that personal taste. As the patient provides all necessary information, I have to see him, have to listen to him, I have to put myself in the patient’s position to find out how individually I am allowed to design the personal work. When creating combination restorations I elaborate aesthetics and set up **together with the patient** at a very early stage. This can be realised very well and easily with the Creopal full contour teeth and Shells. Their design allows sufficient flexibility for individualisation. When I reduce e.g. the incisal area, the tooth will be more chromatic; without much effort I can thus create even more naturalness. There is **a lot of potential** in shape and structure of the Creopal **for individual tooth aesthetics**. It would be a pity taking it from the board without shaping it.

After the set-up I complete the gingiva portions with composite and individualise the teeth with stains – done! To reconstruct the gingiva, I cover the transition between framework and natural gingiva with opaque paste material and layer more and more transparent to the outside. But I never layer too uniformly or too harmonically. Nature is rarely symmetric. Thus I always layer in a certain controlled chaos. My motto: you should be more courageous using different shades, otherwise the effects will disappear very quickly.

In the following **combination prosthesis with Creopal Shells and Creation VC Gingiva** I supplied an edentulous patient with a removable combination prosthesis on four implants. He is a craftsman and a restoration not too bright and symmetric with distinct traces of use, food and beverages fits very well. In this respect I intentionally included a “craftsmanship colour character” in my work – in a controlled, creative way.

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Fig. 3: The framework is conditioned with Metal Primer after blasting and blow-drying (no steam cleaner!).

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
Fig. 4: The Creopal VC Shells are individually adjusted to the framework structure.

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
Fig. 7: The space available between the framework structure and the Creopal Shells is controlled with the aid of the matrix.

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Fig. 8: After sandblasting with aluminium oxide and cleaning of the shells, the Creopal Primer is applied and subsequently sufficiently light-cured.

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Fig. 11: The restoration is polymerised in the light-curing unit; subsequently the excesses are reduced.

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Fig. 12: 1 to 2mm gingiva opaque is applied palatinally at the transition to the metal framework as preparation of the gingiva layering.

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
Fig. 13: An opaque gingiva material is applied to “block” the transition between the natural gingiva and the metal framework.

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Fig. 14: The palatal transition “gingiva/metal” is also covered with G6.

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
Fig. 15: Triangles are applied with the G5 material to darken the interdental areas.

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Fig. 16: Labial over-layering in the apical area from the gingival margin to the natural gingiva.

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

Fig. 17: The flexible parts of the gingiva are individualised with stains and gingiva modifiers – not uniformly, but with a “creative chaos”.

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Fig. 18: The gingiva is individualised also in the palatal area with the Make up Neo stains Illusion and Gingiva as well as the gingiva modifiers.

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Fig. 25: Traces of use, food and beverages are applied with the stains in the cervical, interdental and occlusal area.

- Make up Neo MUN-4 mandarine
- Make up Neo MUN-8 medium brown
- Make up Neo MUN-IL Illusion
- Make up Neo MUN-12 grey

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Fig. 26: The dental prosthesis is polished with customary polishing agents, bison polishing brushes, linen buffs etc.

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Fig. 27: Pay special attention to smooth transitions to the palate anatomy when polishing the palatal areas.

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Fig. 28: Creopal Shells and Creation VC Gingiva – the perfect combination of white and red!

Abb. 29: Die fertige implantatprothetische Kombiarbeit.
Fig. 29: The finished implant-based combination prosthesis.

Vertreiber | Distributor
Creation Willi Geller International GmbH
Koblacherstraße 3, 6812 Meiningen, Austria
Tel. +43 (0)5522 76784
Fax +43 (0)5522 76784-59
info@creation-willigeller.com
www.creation-willigeller.com

Fotografie | Photography
Creation Willi Geller
Nando Aeschlimann
Robert Arvai
Uwe Gehringer
Annette v. Hajmasy
Michael Zangl

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Tooth shades A1 – D4 according to the Vita® shade guide. Vita® is a registered trade mark of VITA Zahnfabrik H. Rauter GmbH & Co. KG, Bad Säckingen, Germany.